

Celemics Whole Exome Sequencing Panel

Key Features

- Complete Whole Exome Coverage
- Superior Performance in the Market
- Gene Add-On Service
- FASTQ to Clinical Interpretation Capability
- Rapid Same-Day Workflow
- No Need for Heavy Instruments
- Complete Walkaway Automation
- Flexible Integration with NGS Sequencers

Celemics Whole Exome Sequencing (WES) Panel is a comprehensive solution that covers all target regions of major WES panels available in the market.

With a target size of 37.1 Mb, the panel does not compromise performance in terms of coverage and uniformity, enabling highly efficient and cost-effective sequencing of the human whole exome. The panel coverage spans across exon regions from RefSeq, CCDS, and GENCODE.

Celemics panels also perform well against hard-to-capture regions such as GC-rich regions. Regarding spike-in options, we have the ability to customize our WES panel according to your specific needs by including mitochondrial or intronic regions upon request.

The panel is also fully supported by Celemics Analysis Service (CAS), our end-to-end bioinformatics solution.

Complete Whole Exome Coverage

What differentiates the Celemics Whole Exome Sequencing Panel from other WES panels in the market?

Most researchers look for complete coverage especially when it comes to WES. However, most WES panels in the market vary in their target regions and some compromise with coverage, even deleting hard-to-capture regions in order to enhance performance. Celemics has developed a WES panel that covers the regions of all four major WES panels in the market, spanning the coding regions from RefSeq, CCDS and GENCODE. The Celemics WES Panel provides the most comprehensive coverage of protein-coding regions, thereby enabling marker discovery to diagnostics.




Figure 1. Complete Whole Exome

The Celemics WES Panel covers the target regions of all major whole exome panels in the market, which include WES panels from company A, company I and company T. Among the regions that Company A and T failed to cover from their own target regions, Celemics covers an additional 60 Kb and 306 Kb of Company A and T respectively (data not shown). Most of these regions are challenging to capture due to GC-rich and repeated sequences. Through Celemics' optimized probe design and panel synthesis technology, the Celemics Whole Exome Panel is able to successfully cover the challenging regions that other company products struggle with.




Figure 2. Superior Capture Performance and Coverage

A mutation in PAFAH1B1 causes Isolated Lissencephaly Sequence (ILS) & Miller-Dieker syndrome. While other competitor panels fail to capture the A-T rich regions in this gene, the Celemics WES Panel successfully covers the region.




Figure 3. Exhaustive Coverage for Each Gene

The Celemics WES Panel covers each gene with thorough coverage in comparison to competitor products. The bar graphs indicate the percentage of genes that are covered at (A) 20X depth and (B) 30X depth. The data from the three panels were downsampled to 5.4 Gb. (C) The IGV figure demonstrates the superior coverage performance of the Celemics Whole Exome Panel against the TIGD1 gene compared to other competitor products.

Exceptional Target Capture Performance

Celemics provides market-leading target capture performance due to probe design and reagent optimization technology. Despite some companies who resort to masking the hard-to-capture regions (such as GC- or AT-rich regions and homologous regions) or completely omit the regions from their target in order to enhance the result quality, Celemics provides both high coverage and on-target ratio without reducing the number of target regions. With Celemics' proprietary technology, the Celemics WES Panel captures regions that no other companies could capture with quality coverage and uniformity. The all-around performance of Celemics' WES panel, regardless of sample type, allows for highly sensitive, cost-effective and time-saving sequencing of the whole exome.




Figure 4. Superior Performance in the Market

Celemics WES Panel shows exceptional performance compared to other competitor products in both germline and FFPE analysis, measured by (A) 1-PCT OFF BAIT (higher the better), (B) Fold 80 base penalty (lower the better). Third-party laboratories (Certified Service Providers) conducted a comparison study between the Celemics WES Panel, Company A and Company I panels. Reference materials NA12878, NA12891, NA12892 for germline analysis, and HD832 for FFPE analysis were used with same amount. Illumina instruments were used for the sequencing. The data from the three panels were downsampled to 5.4 Gb.




Figure 5. Exceptional Uniformity across Low and High GC Regions

(A) The Celemics WES Panel demonstrates minimal deviation, yielding 0.166 standard deviations (lower the better) across GC-rich and AT-rich regions in comparison to competitor products yielding 0.199 and 0.356 standard deviations. (B) The bar graphs shown in different GC ratios also illustrate the consistent uniformity of Celemics WES Panel in comparison to the competitor products.

Simple Gene Add-On Customization




Figure 6. Flexible Customization with Gene Add-On Service

Celeemics provides client-specific customization which is normally not available in the market. Most whole exome panels in the market are provided as a pre-set “ready-to-use” product and options for capturing the intronic or mitochondrial regions are often limited. Fortunately, Celeemics is highly experienced with customized panels, allowing us to offer the Celeemics WES bundled with a customizable Gene Add-On Service option. The addition of intronic regions or mitochondrial genomes to the WES panel enhances comprehensive assessment by capturing clinically significant genes embedded within those regions. The customized WES panel is provided with an option for evaluation through an in-house performance test for best quality results at an additional price.

Rapid Same-Day Workflow

Although hybridization capture has great advantages including minimized bias, stable and reliable data results from a variety of sample types, the complexity of the workflow and the long prep time have been obstacles to the users. Celeemics has developed a new workflow and have incorporated it into our WES package to significantly simplify the process and reduce the experiment time. The conventional method took 2-3 days to complete one sequencing experiment. Now with Celeemics’ newly developed method, the whole experiment and the NGS run can be started on the same day.

Conventional workflow
» 20 hours

Same-day workflow
» 5-8 hours




Figure 7. Newly Developed Same-Day Workflow

The figure demonstrates that Celeemics has significantly reduced the time for performing Whole Exome Sequencing from the conventional 20 hours to 5 hour minimum workflow.

No Need for Heavy Instruments

In order to perform Library Preparation prior to Target Enrichment and Sequencing, it is often required to have heavy instruments (such as a vacuum concentrator or sonicator, etc.), which are barriers against complete automation. Even with an automation protocol, using these heavy instruments is inevitable and is often burdensome to the users. Celemics has successfully eliminated the need for heavy instruments by substituting them with a more convenient solution of enzymes and beads. After rigorous validation that consistently showed reliable performance, we have optimized this workflow to enable the benefit of a complete walkaway solution.




Figure 10. No Need for Heavy Instruments

Celemics has developed and optimized an enzymatic fragmentation, bead-based concentration, and normalization process, which eliminates the requirements for heavy instruments such as a sonicator (mechanical fragmentation), vacuum concentrator, size and amount characterization device, etc.

Flexible Integration with NGS Sequencers & Complete Walkaway Automation




Figure 11. Compatibility with NGS Sequencers and Automation

The Celemics WES Panel is seamlessly integrated with all NGS instruments from Illumina, MGI, and Ion Torrent. Since there are no heavy instruments required, the experiment can be carried out with complete automation.

Ordering Information



| Product Sorting Category / Sub-category | Applied Platform | Package Option | Product Unit |
|---|----------------------------|--|--|
| Ready-to-use panels / Whole Exome Sequencing | Illumina, Ion Torrent, MGI | 1) All-in-one Package 2) Standard Package 3) Target Enrichment Package | 1) Reaction basis: 16, 48, 96 rxn 2) Sample basis: 16, 96, and more options *Pre-capture pooling options: 4, 8, 12 |

*Promotion: pre-capture pooling of 4 sample libraries

Celemics, Inc.

ADD 19F - 20F, Bldg. A, BYC Highcity, 131, Gasandigital 1-ro, Geumcheon-gu, Seoul, 08506, Korea
TEL +82-2-6746-8067 | **FAX** +82-2-6746-8073 | **E-MAIL** support@celemics.com
www.celemics.com

